A Solution to VEPTR Iliac S-Rod Migration

Richard E. McCarthy, M.D.

Professor, Departments of Orthopaedics and Neurosurgery

Arkansas Children's Hospital

University of Arkansas for Medical Sciences

Frances L. McCullough, MNSc, RNP
Arkansas Children's Hospital
University of Arkansas for Medical Sciences

Introduction

- Stable fixation to the pelvis is an essential part of rib to pelvis VEPTR constructs
- A single S-rod is commonly used for pelvic fixation
- Migration of the S-rod into the body of the ilium has been a common problem
- A modification using dual S-rods joined by a domino has been used to provide a more stable anchor point



Methods

- Charts and radiographs were reviewed on all patients who had placement of VEPTR's with rib to pelvis constructs
- A construct of two short S-rods is used over the posterior iliac crest.
 One is directed medial while the other is directed lateral with a domino connector locking them together. A second domino above is attached to the longer of the S-rods for vertical rod attachment





Results


- Seven patients
- Diagnoses:
 - Syndromic 5
 - Neuromuscular 2
- Average age at placement
 6+8 years of age
 - (range 1+11 years to 13+11 years)
- Average length of follow-up
 - 1+2 years
 - (range 6 weeks to 2+10 years)
- Average number of lengthenings per patient 1.5
 - (range 0-6)





Complications

- One broken rod at the pelvis which was revised
- One skin breakdown
 - 2 year old arthrogryposis





Conclusion

 The use of dual S-rods at the posterior iliac crest offers a stable fixation without migration or pull-out at 14 month follow-up